

# Custom Picosecond Amplifiers

## PicoFlux Custom



### APPLICATIONS

- ▶ Multi-stage OPCPA pumping
- ▶ Non-linear optics
- ▶ Other spectroscopic and nonlinear optics applications...

### FEATURES

- ▶ Flash lamp or diode pumped multi-channel PicoFlux systems
- ▶ Each of the channels can be **tailored according to pumping requirements**
- ▶ High energy PicoFlux version with **variable burst**
- ▶ Hybrid PicoFlux with fiber front-end and **Yb:YAG** amplifiers – **1 ps, 8 mJ** output at **10 kHz**
- ▶ **1 – 300 ps** pulse duration
- ▶ From Single Shot to 10 kHz pulse repetition rate
- ▶ Internal or external seeding source
- ▶ Advanced beam shaping for high pulse energy
- ▶ Thermally induced birefringence compensated
- ▶ Low jitter synchronization pulses below 10 ps RMS jitter
- ▶ Vacuum image relay system
- ▶ Optional temperature stabilized second, third and fourth harmonic generators

Multiple Channel PicoFlux series picosecond amplifiers were designed and manufactured for multiple stage OPCPA pumping. Systems can be specially tailored for customer's needs and have up to 8 pumping channels with different wavelength, energy, pulse duration, spatial and temporal profiles, adjustable delay, image translation to customers specified location and various other features. Short pulse duration, excellent pulse-to-pulse stability, superior beam quality makes PicoFlux series picosecond amplifiers well suited for other applications as well.

### Regenerative amplifier / Power amplifier design

PicoFlux series amplifiers are designed to be seeded by external seeding source. Diode pumped regenerative amplifier ensures amplification of seed signal to stable mJ level pulse for amplification in linear amplifiers. Advanced beam shaping ensures smooth, without hot spots beam spatial profile at the laser output. Low light depolarization level allows high efficiency generation of up to 4th harmonic with optional build-in harmonic generators.

Alternatively EKSPLA can offer an internal seeder meeting customer's requirements.

### Full-fiber front-end

Novel EKSPLA developed fiber front-end opens up new set of unique features for PicoFlux systems like dual wave seeder that also has burst formation functionality due to active fiber loop technology. The front-end also permits to offer powerful, ultrafast and higher frequency Ytterbium amplifiers for wider opportunities in OPCPA pumping and other scientific uses.

## SPECIFICATIONS

| Model  | P2k10-x4  | P1301k-x8   | P1k100-Burst                              | P810k-1030                            |
|--|---|---|---|---------------------------------------|
| <b>MAIN SPECIFICATIONS <sup>1)</sup></b>       |   |   |   |                                       |
| Output energy                                  |   |   |   |                                       |
| Fundamental                                    | 4 × 2 200 mJ <sup>2) 3)</sup>                             | 8 × 130 mJ  | 1.4 J Burst<br>(4×300 mJ + 4×50 mJ)       | 8 mJ @ 1030 nm                        |
| SH output <sup>4) 5)</sup>                     | 4 × 1 400 mJ  | 8 × 85 mJ   | NA  | NA                                    |
| TH output <sup>4)</sup>                        | 4 × 660 mJ  | 8 × 50 mJ   | NA  | NA                                    |
| FH output <sup>4)</sup>                        | 4 × 220 mJ  | 8 × 15 mJ   | NA  | NA                                    |
| Pulse repetition rate                          | 10 Hz   | 1 kHz   | 100 Hz                                    | 10 kHz                                |
| Pulse duration <sup>6)</sup>                   | 90 ± 10 ps  | 90 ± 10 ps  | 90 ± 10 ps                                | 1 ± 0.2 ps                            |
| Pulse energy stability <sup>7)</sup>           |   |   |   |                                       |
| Fundamental                                    | ≤ 0.6 %   | ≤ 0.5 %   | ≤ 1 %                                     | ≤ 0.5 %                               |
| SH output <sup>4)</sup>                        | ≤ 0.8 %   | ≤ 0.8 %   | NA  | NA                                    |
| TH output <sup>4)</sup>                        | ≤ 2 %   | ≤ 2 %   | NA  | NA                                    |
| FH output <sup>4)</sup>                        | ≤ 3 %   | ≤ 3 %   | NA  | NA                                    |
| Long-term power drift <sup>8)</sup>            | ± 2 %   | ± 1.5 %   |   |                                       |
| Beam spatial profile                           | Super-Gaussian <sup>9)</sup>                              | Super-Gaussian <sup>9)</sup>                              | Super-Gaussian and Gaussian <sup>9)</sup> | Gaussian <sup>9)</sup>                |
| Beam diameter <sup>10)</sup>                   | ~ 23 mm   | ~ 7 mm  | ~ 11 & 5 mm                               | ~ 6 mm                                |
| Beam pointing stability <sup>11)</sup>         | ≤ 30 μrad   | ≤ 20 μrad   | ≤ 20 μrad                                 | ≤ 20 μrad                             |
| Beam divergence                                | ≤ 0.5 mrad  | ≤ 0.5 mrad  | ≤ 0.5 mrad                                | ≤ 0.5 mrad                            |
| Pre-pulse contrast <sup>12)</sup>              | > 200:1   | > 200:1   | > 200:1                                   | > 200:1                               |
| Optical pulse jitter <sup>13)</sup>            |   |   |   |                                       |
| Trig out                                       | ≤ 100 ps  | ≤ 100 ps  | ≤ 100 ps                                  | ≤ 50 ps                               |
| Pre-Trig out                                   | ≤ 50 ps   | ≤ 50 ps   | ≤ 50 ps                                   | ≤ 50 ps                               |
| With -PLL option                               | ≤ 2 ps  | ≤ 2 ps  | ≤ 2 ps                                    | NA                                    |
| Polarization                                   | Linear  | Linear  | Linear                                    | Linear                                |
| <b>PHYSICAL CHARACTERISTICS <sup>14)</sup></b> |   |   |   |                                       |
| Laser head size (W×L×H mm)                     | 1500 × 3600 × 500, 2 pc.                                  | 1500 × 3600 × 500, 4 pc.                                  | 700 × 2000 × 300                          | 900 × 1200 × 300                      |
| Power supply size (W×L×H mm)                   | 553 × 600 × 1800, 4 pc.                                   | 553 × 600 × 1800, 4 pc.                                   | 553 × 952 × 600                           | 553 × 952 × 600                       |
| Umbilical length <sup>15)</sup>                | 5 m   | 2.5 m   | 2.5 m                                     | 3 m                                   |
| <b>OPERATING REQUIREMENTS <sup>16)</sup></b>   |   |   |   |                                       |
| Electrical power                               | 208, 380 or 400 VAC, three-phase, 50/60 Hz <sup>17)</sup> | 208, 380 or 400 VAC, three-phase, 50/60 Hz <sup>17)</sup> | 208 – 240 VAC, single-phase, 50/60 Hz     | 208 – 240 VAC, single-phase, 50/60 Hz |
| Power consumption <sup>18)</sup>               | ≤ 40 kVA  | ≤ 60 kW   | ≤ 5 kW                                    | ≤ 3.5 kW                              |
| Water supply                                   | ≤ 40 l/min, 2 Bar, max 15 °C                              | ≤ 40 l/min, 2 Bar, max 15 °C                              | ≤ 5 l/min, 2 Bar, max 15 °C               | ≤ 15 l/min, 2 Bar, max 15 °C          |
| Operating ambient temperature                  | 22 ± 2 °C   | 22 ± 2 °C   | 22 ± 2 °C                                 | 22 ± 2 °C                             |
| Storage ambient temperature                    | 15 – 35 °C  | 15 – 35 °C  | 15 – 35 °C                                | 15 – 35 °C                            |
| Relative humidity (non-condensing)             | ≤ 80 %  | ≤ 80 %  | ≤ 80 %                                    | ≤ 80 %                                |
| Cleanliness of the room                        | ISO Class 7   | ISO Class 7   | ISO Class 7                               | ISO Class 7                           |

<sup>1)</sup> Due to continuous improvement, all specifications are subject to change without notice. The parameters marked 'typical' are indications of typical performance and will vary with each unit we manufacture. Presented parameters can be customized to meet customer's requirements. All parameters measured at 1064 nm if not stated otherwise.

<sup>2)</sup> 2 200 mJ energy is achieved with Super-Gaussian spatial beam profile of 11<sup>th</sup> or higher order (with steep edges). If lower order Super-Gaussian is required maximum pulse energy will be limited to 2 000 mJ.

<sup>3)</sup> 2 500 mJ output energy is available upon request with longer pulse duration.

<sup>4)</sup> Harmonic outputs are optional. Specifications valid with respective harmonic module purchased. Outputs are not simultaneous.

<sup>5)</sup> Second harmonic specification is valid when only SH option is ordered. If TH/FH options are orders second harmonic efficiency is reduced to ~50 %.

<sup>6)</sup> Standard pulse duration is 90 ps. Other pulse durations can be ordered within range of 20 ps – 300 ps. Shortening the pulse duration below 90 ps will reduce the output energy proportionally.

<sup>7)</sup> Under stable environmental conditions, normalized to average pulse energy (RMS, averaged from 60 s).



- <sup>8)</sup> Measured over 8 hours period after 30 min warm-up when ambient temperature variation is less than  $\pm 2$  °C.
- <sup>9)</sup> Super-Gaussian spatial mode of 6-11<sup>th</sup> order in near field.
- <sup>10)</sup> Beam diameter is measured at signal output at 1/e<sup>2</sup> level for Gaussian beams and FWHM level for Super-Gaussian beams.
- <sup>11)</sup> Beam pointing stability is evaluated as movement of the beam centroid in the focal plane of a focusing element (RMS, averaged from 60 s).
- <sup>12)</sup> 1000:1 contrast available upon request.

- <sup>13)</sup> Optical pulse jitter with respect to electrical outputs:  
– Trig out > 3.5 V @ 50  $\Omega$   
– Pre-Trig out > 1 V @ 50  $\Omega$   
– PLL option > 1 V @ 50  $\Omega$
- <sup>14)</sup> System sizes are preliminary and depend on customer lab layout and additional options purchased.
- <sup>15)</sup> Longer umbilical with up to 10 m for flash lamp pumped and up to 5 m for diode pumped systems available upon request.
- <sup>16)</sup> The laser and auxiliary units must be settled in such a place void of dust and aerosols. It is advisable to operate the laser in air

conditioned room, provided that the laser is placed at a distance from air conditioning outlets. The laser should be positioned on a solid worktable. Access from one side should be ensured.

- <sup>17)</sup> Voltage fluctuations allowed are +10 % / -15 % from nominal value.
- <sup>18)</sup> Required current rating can be calculated by dividing power rating by mains voltage. Power rating is given in apparent power (kVA) for systems with flash lamp power supplies and in real power (kW) for systems without flash lamp power supplies where reactive power is neglectable.

## OPTIONS

| Option     | Description   | Comment   |
|------------|---|---|
| -P20...300 | Custom pulse duration between 20 ps and 300 ps                | Available with internal and external seeder. Shortening the pulse duration below 90 ps will reduce the output energy proportionally   |
| -50/100    | 50 Hz or 100 Hz pulse repetition rate                         | Energy can be increased ~4 times compared to 1 kHz systems  |
| -2k        | 2 kHz pulse repetition rate                                   | Reduces the output energy of fundamental by ~50 %   |
| -G         | Gaussian like spatial beam profile                            | Reduces the output energy of fundamental by ~80 %   |
| -FS        | External seeder input via motorized spectral broadening stage | Requires > 1.5 nJ per pulse @ 800 nm, 100 fs  |
| -PLL       | Phase Lock Loop option for precise lock to external RF signal | Electrical to optical signal jitter $\leq 3$ ps   |
| -SH/TH/FH  | Second, third and fourth harmonic outputs                     | Conversion efficiency from fundamental respectively ~50 %, ~30 % and ~10 %. Harmonic outputs not simultaneous with fundamental output |
| -AW        | Water-to-Air cooling  | Replaces or supplements Water-to-Water cooling unit. Heat dissipation equals total power consumption                                  |

## PERFORMANCE




Fig 1. Typical PicoFlux Ytterbium system amplifier system near field beam profile at 1030 nm (imaged from laser output)




Fig 2. Typical High Energy PicoFlux amplifier system near field beam profile at 1064 nm (imaged from laser output)




Fig 3. Typical High repetition rate PicoFlux amplifier system far field beam profile at 532 nm (imaged from SH crystal)




Fig 4. Typical output pulse spectrum of the PicoFlux Ytterbium system




Fig 5. Typical pulse shape (FROG measurement) PicoFlux Ytterbium system

## STABILITY




Fig 6. Typical long-term energy stability of High repetition rate PicoFlux system




Fig 7. Typical long-term energy stability of High Energy PicoFlux system

## ORDERING INFORMATION

**Note:** Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer than 1 hour then laser (system) needs warm up for a few hours before switching on.

### PicoFlux P(1)(2)-x(3)-(4)

Energy level:  
130 → 130 mJ  
2k → 2200 mJ

Pulse repetition rate:  
10 → 10 Hz  
1k → 1 kHz

Any additional options:  
See 'Options' table

Number of channels:  
x4 → four channels  
x8 → eight channels